


# *LeadCross*<sup>TM</sup> Guide Wire

**LeoMed**  
乐奥医疗 为生命健康服务

---

GI Endoscopic Accessories Series 

---

# LeadCross™ Guide Wire

Used for guiding the catheter or other devices to the digestive tract for interventional diagnostic and surgery endoscopically.

- **Soft tip**  
Easy to insert and less tissue damage.

- **Good Visibility**  
Colorful coating is easy to identify. Convenient for following the device movement under endoscope.

- Unique super-elastic nitinol core designed to promote excellent kink resistance, Dynamic nitinol core intended to promote pushability for easy advancement.

- Hydrophilic coating designed to facilitate crossing of tight lesions.

## Product Specification

REF	Diameter ( inch )	Working length ( mm )	Tip shape	Hardness
L023520SS00	0.035	2000	Straight	Super stiff
L023526SS00	0.035	2600		
L023545SS00	0.035	4500		
L023520SN00	0.035	2000	Straight	Normal stiff
L023526SN00	0.035	2600		
L023545SN00	0.035	4500		